



# BIOSOLVE-IV FIRST COHORT

## 36-month follow-up of **first cohort** with **1,075 patients**<sup>1</sup>

### Conclusions

- **Magmaris**<sup>®</sup> Resorbable Magnesium Scaffold (RMS) shows a good safety profile of the scaffold in a real-world setting up to 36 months for the first cohort.<sup>1</sup>
- The Target Lesion Failure (TLF)\* rate of 8.2% at 36 months for the first cohort is comparable to low TLF rates of contemporary drug-eluting stents (DES)<sup>2,3</sup> over the same time period.
- Only 1 very late scaffold thrombosis (ST)\*\* case occurred up to 36 months for the first cohort of BIOSOLVE-IV.

### Study design

Prospective, multi-center, real-world setting registry

### Patients

19.2% of NSTEMI patients included


### Endpoints

#### Primary endpoint

- TLF\* at 12 months

#### Secondary endpoints

- Clinically-driven Target Lesion Revascularization (Clinically-driven TLR) and Target Vessel Revascularization (TVR)
- Cardiac death
- Target-Vessel Myocardial Infarction (TV-MI)
- Definite or probable ST at 12 months powered to show the superiority of Magmaris vs. historical data of Absorb for full cohort based on a one one-sided exact binomial test
- Procedure and device success



### Patient characteristics

	n = 1,075	
Age, yrs <sup>o</sup>	61.3	± 10.5
Male	806	75.0%
Hypertension	724	67.3%
Hypertlipidemia	713	66.3%
Smoking	654	60.8%
Diabetes mellitus	228	21.2%
Insulin dependent	45	19.7%
Non-insulin dependent	183	80.3%
History of MI	219	20.4%
Previous percutaneous coronary intervention	287	26.7%
<b>NSTEMI</b>	<b>206</b>	<b>19.2%</b>

### Lesion location

	n	
LAD	561	50.0%
LCx	213	19.0%
RCA	333	29.7%
Ramus intermedius	14	1.2%

### Lesion characteristics


	n	
Lesion length (mm) <sup>o</sup>	14.9	± 4.2
Reference vessel diameter (mm) <sup>o</sup>	3.2	± 0.3
<b>AHA/ ACC lesion class B2/C</b>	<b>170</b>	<b>15.2%</b>
Calcification moderate/severe	82	7.3%
<b>Bifurcation lesions</b>	<b>57</b>	<b>5.1%</b>

\* TLF is defined as a composite of Cardiac Death, TV-MI, emergent Coronary Artery Bypass Grafting (eCABG), and Clinically-driven TLR. Peri-procedural MI according to SCAI definition and spontaneous MI according to the extended historical definition.

\*\* At 1,043 days post procedure: Patient presented with inferior STEMI. Thrombotic occlusion of the previously scaffolded segment in the proximal RCA (about 3 years before). TLR with a DES performed.

<sup>o</sup> Data shown as mean ± SD

## TLF\* of Magmaris<sup>1</sup> up to 36 months



## TLF of contemporary DES<sup>2,3</sup> up to 36 months

BIORESORT <sup>2</sup> TLF			BIOFLOW-V <sup>3</sup> TLF		
SES	EES	ZES	SES	DP EES	
n = 1,169	n = 1,172	n = 1,173	n = 884	n = 450	
77	6.7%	86	7.5%	96	8.3%
			70	8.2%	59
					13.6%

The event rates, expressed as n (%), were calculated with the use of the Kaplan-Meier method at 36 months.

DP EES = Durable-polymer everolimus-eluting stent  
 EES = Everolimus-eluting stent  
 SES = Sirolimus-eluting stent  
 ZES = Zotarolimus-eluting stent


## BIOSOLVE-IV First cohort<sup>1</sup> TLF\* and selected components

	BIOSOLVE-IV <sup>4</sup> 12-month FUP n = 1,071 <sup>o</sup>	BIOSOLVE-IV <sup>5</sup> 24-month FUP n = 1,071 <sup>o</sup>	BIOSOLVE-IV <sup>1</sup> 36-month FUP n = 1,071 <sup>o</sup>			
<b>TLF*</b>	<b>45</b>	<b>4.3%</b>	<b>70</b>	<b>6.6%</b>	<b>86</b>	<b>8.2%</b>
Cardiac death	2	0.2%	5	0.5%	10	1.0%
TV-MI	12	1.1%	16	1.5%	21	2.0%
Clinically-driven TLR	41	3.9%	62	5.9%	72	6.8%
Emergent CABG	0	0.0%	0	0.0%	0	0.0%
<b>Scaffold Thrombosis (definite or probable)</b>	<b>5</b>	<b>0.5%<sup>Δ</sup></b>	<b>5</b>	<b>0.5%<sup>Δ</sup></b>	<b>6</b>	<b>0.6%<sup>§</sup></b>

n-values are based on BIOTRONIK data on file and assessed by Kaplan-Meier failure estimate analysis.

<sup>o</sup>Rational for the n-value calculation: 1,075-4 devices not implanted = 1,071

## Def./ prob. ST up to 36 months



BIOSOLVE-IV first cohort demonstrates consistently low scaffold thrombosis rates up to 36 months.<sup>1</sup>

<sup>§</sup> 0.5% scaffold thrombosis rate excluding cases with early antiplatelet or anticoagulant interruption.

<sup>Δ</sup> 4 out of 5 (0.1%) scaffold thrombosis cases had early antiplatelet or anticoagulant interruption after procedure

## Coordinating investigators

Stefan Verheye, MD, PhD, Antwerp Cardiovascular Center, ZNA, Middelheim, Antwerp, Belgium and  
 Michael Kang-Yin Lee, MD, PhD, Queen Elizabeth Hospital, Hong Kong

## Principal investigator

Torzewski Jan, Prof, Cardiovascular Center Oberallgaeu-Kempton, Oberallgaeu, Germany

1. Torzewski J. Safety and performance of Magmaris at 36-months: BIOSOLVE-IV first cohort. Presented at: EuroPCR; 2022; ClinicalTrials.gov: NCT02817802; 2. Buiten A et al. Thin, Very Thin, or Ultrathin Strut Biodegradable - or Durable-Polymer-Coated Drug-Eluting Stents 3-Year Outcomes of BIO-RESORT. JACC: CARDIOVASCULAR INTERVENTIONS. 2019. doi: 10.1016/j.jcin.2019.04.054; 3. Kandzari D.E et al. Ultrathin Biodesorbable-Polymer Sirolimus-Eluting Stents Versus Thin Durable-Polymer Everolimus-Eluting Stents for Coronary Revascularization 3-Year Outcomes From the Randomized BIOFLOW V Trial. JACC: CARDIOVASCULAR INTERVENTIONS. 2020. 13 (11): 1343-53. doi:10.1016/j.jcin.2020.02.019; 4. Verheye S, Włodarczak A, Montorsi P, et al. BIOSOLVE-IV-registry: Safety and performance of the Magmaris scaffold: 12-month outcomes of the first cohort of 1,075 patients. Catheter Cardiovasc Interv. 2020; 1-8. doi.org/10.1002/ccd.29260; 5. Torzewski J. Safety and performance of Magmaris at 24 month follow up of BIOSOLVE IV. Presented at: eEuroPCR; 2021; virtual congress. ClinicalTrials.gov: NCT02817802.

All events have been adjudicated by an independent clinical event committee. BIOSOLVE-IV registry is based on Kaplan-Meier failure estimate analysis including censored observations. Please consult IFU for indications and use in special populations.

Magmaris is a trademark or registered trademark of the BIOTRONIK Group of Companies.

BIOTRONIK AG  
 Ackerstrasse 6  
 8180 Bülach  
 Switzerland

Tel +41 (0) 44 8645111  
 Fax +41 (0) 44 8645005  
 info.vi@biotronik.com  
 www.biotronik.com

© 2022 BIOTRONIK AG – All rights reserved.  
 Specifications are subject to modification, revision and improvement.

**BIOTRONIK**  
 excellence for life